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Thank you for purchasing a Casada product.

These instructions for use contain important safety and care instructions and explain how to use the product. Please read them carefully and keep the instructions in a suitable location for future reference.

SECURITY INSTRUCTIONS

- Keep the Reflex Disc away from heat and sun.
- Wipe the Reflex Disc regularly using a dust-free, lightly moistened cloth.
- Avoid contact with sharp objects such as rings, sharp fingernails or others.
- An improper use of the band will result in a loss of the liability coverage.

WARNING: A plastic bag was used to pack this product. Keep this plastic bag away from babies and children to avoid suffocation. Plastic bags are no toys.

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Please step on the Reflex Disc; you can move your feet on the mat randomly.

Suggestion: 10 minutes each time. This action can be repeated many times.

Benefit: Improves blood circulation and promotes "Chi". Regular use promotes the flow of energy.

Suitable for: Elderly people, people with a lack of exercise, cold feet and muscle discomfort.

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Please step on the Reflex Disc using your toes, you can move on the mat to the left or to the right.

Suggestion: 5 minutes each time. This action can be repeated many times.

Benefits: Relieves shoulder- and headaches, helps to relieve stress.

Suitable for: Housewives, persons with sleep disorders, stressed adults.

USAGE

Place the Reflex Disc on an evenly and flat surface.

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Please step on the Reflex Disc using your heels and move on the mat to the left or to the right.

Suggestion: 10 minutes each time. This action can be repeated many times.

Benefits: Good for digestion system.

Suitable for: Middle-aged persons, working adults, people who are "Chi" – less.

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Use the Reflex Disc while working, watching TV or reading. Place both or just one foot gently and enjoy the soothing effect.

Suggestion: 20 minutes or more. This action can be repeated many times.

Benefits: Improves blood circulation and relieves digestive problems.

Suitable for: Working adults in air-conditioned environment, elderly people and those with overweight and lack of exercise.

TECHNISCHE DATEN TECHNICAL DATA

Maße measure:

Ø 33,5 cm

Gewicht weight:

0,45 kg

Material material:

EVA-Kunststoff

Zertifikate certificates:

CASADA
healthcare